



## ECHELON Smart Plus

THE ENLIGHTENED CHOICE FOR HIGH  
PRODUCTIVITY AND DIAGNOSTIC POWER



the next level in magnetic resonance imaging



## KEY FEATURES


- SmartCOMFORT
- SmartSPEED
- SmartECO
- SmartSPACE

# **ECHELON Smart Plus heralds the dawn of a new standard for 1.5T superconductive MRI**

The ECHELON Smart Plus features a small footprint with economics that do not compromise diagnostic quality and speed.

Based on FUJIFILM's proprietary technology, this system has opened up the potential for 1.5T systems, providing superb image quality and superior install flexibility inherited from permanent-magnet MRI systems.

ECHELON Smart Plus offers new options for superconductive MRI.



**96%** REDUCTION\*

# SmartCOMFORT

## Patient-friendly quiet examinations

### Quiet examinations with SoftSound

Various technologies exist to reduce MRI acoustic noise. However, reducing acoustic noise often compromises image quality and may extend scan time, making it unsuitable for routine examination. Other approaches need dedicated hardware. FUJIFILM's SoftSound noise reduction technology reduces the acoustic noise by up to 96%\* maintaining image quality and scan time, and without any dedicated hardware.

\* Compared to ECHELON Smart without SmartCOMFORT and in dependence of the sequence parameters

#### Plus

IP-RAPID can be used along with SoftSound



SoftSound + IP-RAPID

#### Plus

SoftSound can also be used with the "BSI" high-speed 3DT2\* weighted images and "3D GEIR"/"3D RSSG" high-speed 3DT1 weighted images.




BSI

3D GEIR

3D RSSG

# Experience enhanced SPEED with FUJIFILM's SynergyDrive


## Conventional equipment



Conventional equipment

Enter MRI room

Patient setting



## Workflow coil system

Positioning, imaging and post-processing

ECHELON Smart Plus

Enter MRI room

Patient setting

Patient positioning

### ■ AutoPose

The new Scanogram function allows information to be displayed rapidly on cross-sectional areas of regions for operators.

FUJIFILM's SynergyDrive offers enhanced speed during MRI examinations without reducing image quality or increasing acoustic noise.

### FUJIFILM's SynergyDrive offers efficient workflow

"SynergyDrive" offers various functions and applications to simplify the series of operation procedures from entering and leaving the examination room to shortening examination time. This is also expected to improve management efficiency.

## ECHELON Smart Plus

**Patient positioning****Scanning**

### AutoExam

are all completed in one step. The examination time is shortened by simplifying operations.

**Scanning****Post-processing****Exit MRI room****■ IP-RAPID****■ AutoClip**

Digital image processing function that automatically creates clipping images of head MRA images.

**■ DICOM transfer**

Automated image transfer by DICOM

**■ High speed Prescan**

Prescan, which is performed as a preparatory procedure, is also enhanced by reducing the time before the main scan.



### IP-RAPID A cutting-edge high-speed imaging technology

IP-RAPID is high-speed imaging using iterative processing, reducing scan times by up to 60% while maintaining image quality comparable to conventional imaging methods.

**Up to****60 % REDUCTION\***

\* Compared to ECHELON Smart without SynergyDrive




# Experience newly enhanced

## SPEED with FUJIFILM's IP-RAPID


Normal scan

Total 16:42



With IP-RAPID

Total 8:14



### IP-RAPID enables both high image quality and reduced scan time

IP-RAPID is a new technology that can reduce scan time while maintaining image quality.

By optimizing and combining undersampling and iterative reconstruction, it can be combined with various regions and functions.

SNR and spatial resolution can also be improved while maintaining the scan time.

Total scan time



**50%**  
REDUCTION\*

\* Value is compared between the following examples, base on the left is ECHELON Smart without IP-RAPID

## SynergyDrive

Workflow coil system    AutoExam  
■ AutoPose

Enter MRI room

Patient setting

Patient positioning

Scanning

■ AutoClip ■ DICOM transfer


Post-processing

Exit MRI room


IP-RAPID

**IP-RAPID covers a variety of pulse sequences for all anatomical regions.**


Emergency imaging Total 3:59




DWI Total 4:00




Multi-contrast imaging Total 5:54




Gated imaging (ECG, respiratory)



Combined with SoftSound



High-resolution imaging





# Experience enhanced SPEED

## with FUJIFILM's AutoExam



### Auto operation



#### ■ AutoPose brain

AutoPose is a slice line setting support function. When the Scanogram imaging ends, the slices are automatically positioned and allow to register up to five types of presets. Below you see two different examples.

#### ■ Examples of preset screens




### Fully automatic function of MRI examination with simplified operations

AutoExam enables setting of imaging conditions, positioning, image processing, image display, image storage, and image transfer functions to be carried out in one step during examination.


The operator can choose to perform semi-automated examinations and may stop, correct and restart imaging.

## SynergyDrive



### ■ AutoClip

When AutoExam is executed, clipping is performed automatically after MRA imaging.  
It is possible to perform additional clipping on images after automatic clipping.



# SmartECO

Ecological and  
economical running costs

17% LESS  
RUNNING  
COSTS\*

\* Compared to ECHELON Smart




**FUJIFILM's SmartECO effectively reduces power consumption**


A compact scan room with an MRI unit with less heat discharge also reduces the air conditioning requirements of the examination and equipment rooms. The energy saving function together with the reduction of heat discharge can reduce running costs by 17%\*.

\*Depending on operating conditions and other factors.

**without SmartECO**



**with SmartECO**




### Minimize unnecessary power consumption

Superconductive MRI systems generally command high running costs. These costs are mainly related to the high power consumption of the cooling system necessary to maintain superconductivity. ECHELON Smart Plus is equipped with an energy saving function that can stop the cooling system for a certain length of time during periods of non-use or on non-consultation days. This function effectively reduces the power consumption whilst maintaining zero helium boil-off. Furthermore, as the heat emission from the cooling system itself also decreases during these periods, the power consumption of its heat-dissipating unit is also cut.

# SmartSPACE

## Small footprint flexible layout

ECHELON Smart Plus




### Flexible layout even in limited space

It is often a matter of concern whether there is sufficient space in the equipment room for a superconductive MRI system installation.


ECHELON Smart Plus can remove such obstacles during installation.

**41%** SPACE SAVING\*

ECHELON Smart Plus saves 41% of the installation space in the equipment room



**Example of installation**  
**Example of equipment room layout**  
An equipment room can be built using an open plan





# ECHELON Smart Plus CAPABILITIES

**Attain high image speed in clinical practice  
with the combination of FUJIFILM's applications,  
operations, and hardware.**

## SmartAPPLICATION

All Around RADAR

Plaque imaging

FatSep

AutoPose spine

BeamSat TOF

## SmartIMAGE GALLERY

## SmartHARDWARE

Work coil system

High performance

RF system

16ch receiver system

High performance

gradient system

## Sentinel Analytics



# SmartAPPLICATION

## Efficient applications for enhancing MR imaging

### All around RADAR

#### RADAR is applicable for routine head examinations


RADAR reduces motion artefacts and increases the ease of use by making it compatible with multiple sequences, all receiver coils and arbitrary cross sections. RADAR can be used in combination with high-speed imaging as well. ECHELON Smart Plus's "All Around RADAR" is compatible with TOF sequences, GrE sequences and most of the sequences required for routine brain examinations.



#### Effects of RADAR on TOF MRA and GrE T2\*WI

RADAR has been applied to GrE sequences using a high-precision signal correction technology.


This has enabled the combined use with RADAR for all sequences required in routine brain examinations.



### Plaque imaging

#### Diagnosis of plaque characterization

Diagnosis of carotid artery plaque characterization requires a high T1 contrast MRI image.



The non-gated RADAR-SE method (also known as Radial Scan), maintains a constant TR without being affected by pulsation and can conduct scanning with a high T1 contrast appropriate for diagnosis of plaque characteristics.

SIR Map shows the colour map depending on the signal strength ratio after normalizing the ROI signal strength using reference signal strength.

## FatSep

**FatSep is a type of suppression method and is less influenced by changes in magnetic susceptibility**


Multiple images can be obtained with one scan, using the difference in resonance frequency (chemical shifts of water protons and fat protons).



## AutoPose spine

**Assistance with imaging plane setting allows reduced in operation time**

AutoPose spine is a support function for quick and accurate slice setting. AutoPose processing is executed at the end of the Scanogram, and the scanning positions of the AX / SAG / COR cross-sections are calculated simultaneously. The scanning position of the AX slices are set automatically parallel to the intervertebral disc.





## BeamSat TOF

**Improving visibility in haemodynamic changes**

**Selective MRA - Addition of haemodynamic information added to TOF**


Pencil-beam type pre-saturation (BeamSat) pulses based on the application of local excitation are used in TOF imaging to selectively suppress some of the blood flow signals required for identification of the haemodynamics.

Selective MRA scans with BeamSat pulses, which targets a specified blood flow allow signals to be suppressed and the dominant regions to be clearly identified. BeamSat pulses can be set to arbitrary positions and angles using a dedicated GUI.




# SmartIMAGE GALLERY


## Neuro Vascular




## Body



**Spine**



**MSK**



# SmartHARDWARE

FUJIFILM's hardware technologies for enhancing image quality

## Workflow coil system

Receiver coils support ease of patient setting and offer superb image quality

The number of receiver coils that must be set prior to the examination is minimized to reduce replacement time and effort. With a system designed for ease of use and with the adoption of special receiver coils for individual regions, significant reduction in examination time is attained whilst maintaining high image quality.

### Workflow coil system

for head & neck



for spine



for abdomen



for joints



### Workflow coil system Setting examples



**Spine examination**  
- Head and neck coil (posterior side) +  
WIT spine coil



**Head and neck examination**  
- Head and neck coil +  
WIT spine coil



**Abdomen examination**  
- Head and neck coil [posterior side] +  
WIT spine coil + Flex body upper coil



**Extremity examination**  
- Head and neck coil [posterior side] +  
WIT spine coil + GP flex coil


## High performance RF system

### Powerful RF output contributes to stable RF transmission

ECHELON Smart Plus is equipped with an RF power output of 18 kW. This is sufficient to provide clear images without deterioration of image quality even in the FSE sequence that applies refocus pulses continuously.

## 16ch receiver system

### Increasing receiver coil channels improves sensitivity and stability



## Major hardware specifications

ECHELON Smart Plus is equipped with a high performance hardware system which is essential for creating high quality images.

### High performance hardware

#### Gradient system

Maximum gradient strength    33mT/m  
Maximum slew rate                130T/m/s

#### RF system

Maximum RF output    18kW

#### High Order Shimming System


HOSS

# Sentinel Analytics<sup>\*1</sup>

## Improving the uptime through failure sign diagnosis

Achieving higher uninterrupted system availability and optimizing maintenance costs remain challenges for conventional remote support services for medical devices. FUJIFILM has accumulated and analyzed big data to develop a new system that utilizes its "Failure sign diagnosis service" to launch "Sentinel Analytics", a failure sign diagnostic service for superconductive MRI systems.

With the failure sign diagnosis based on IoT<sup>\*2</sup>, the inspection and parts replacement cycles can be optimized and the system availability can be improved.



### Major features and advantages

#### Constant system monitoring

The Sentinel server monitors the system state 24 hours a day.

#### Automatic notification feature

When the Sentinel server detects either a malfunction or a lowered performance of the system, an alert is automatically reported to the FUJIFILM service site. This helps prevent the occurrence of a malfunction. Furthermore, a corrective measure is quickly taken in the case of a malfunction.


#### Direct connection feature

This feature provides service via direct connection of the service site and your system. To track down the causes of a malfunction, we check artefacts and abnormal images, check image data before reconstruction (raw data) and run test programs on the system.


#### Data security

Features such as encryption of communication data and communication based on mutual authentication are available to protect patient information. Furthermore, the specification does not allow recognition of personal information included in patient lists and images (such as a patient's name, sex, weight, age, and date of birth) on the Sentinel server and the service site.

## Highly accurate failure sign diagnostics



## Conceptual diagram of Sentinel Analytics



\*1 Service contract is required.

\*2 IoT (Internet of Things) : A system in which various devices with communication functions exchange information via the Internet to realize identification, monitoring, and control of such devices.

\*3 M2M (Machine-to-Machine) : A system of direct exchange of information between machines via a network without human intervention.



## ECHELON Smart Plus

- ECHELON, Sentinel, Sentinel Analytics, VASC and FatSep are registered trademarks or trademarks of FUJIFILM Healthcare Corporation in Japan and other countries.
- This brochure may contain descriptions of optional functions and products.
- Specifications and physical appearance may change without prior notice.
- Please refer to the operation manual and the related documents for appropriate use of this product.

**FUJIFILM**

Manufacturer

**FUJIFILM Healthcare Corporation**

2-1 Shintoyofuta, Kashiwa-shi, Chiba, 277-0804, Japan  
[www.fujifilm.com/fhc/en](http://www.fujifilm.com/fhc/en)

Importer and distributor for Europe

**FUJIFILM Healthcare Europe Holding AG**

Sumpfstrasse 13, 6312 Steinhausen, Switzerland  
[www.fujifilm.com/hce](http://www.fujifilm.com/hce)

ECHELON Smart Plus-FF/EU-Version/EN,11/2021/v2/NIK